## HIGHT OF THE LIVING LIGHT REAPOUT UPPATE

Starring:

Will Foreman, Andrzej Szelc, David Schmitz

LArIAT Meeting October 31, 2013


## Highlights since last update...

- Plotted gain from 1p.e. response with PMTs using LED pulser + WArP code
- Tried to suppress ringing/noise due to inductance from driving the pulser
  - See "Light readout update (July 18 2013)" on DocDB
- Mounted a SiPM to circuit and began testing the output


### **SER from PMTs**


DAQ/processing code gives reliable single photoelectron response


Increases with HV as expected



## Mounting the SiPM











### Testing the SiPM



Cannot see photoelectron peaks in the SiPM like we can in the PMTs...

 Tried amplifying using NIM module for positive PMT amplification (Phillips 774)
 → amps noise as well as signal

### Testing the SiPM



Bright and quick LED flash...


- Seen by PMT as sharp 500mV pulse
- Seen by SiPM as unstable ~10mV pulse with bizarre ringing and long tail

Need better response in order to calibrate the SiPM and be sensitive to SER...

### Next steps

#### Try different SiPM circuit

- forward vs. reverse voltage bias?
- Can we amplify signal without amplifying the noise with it?
- Mount SiPM to an array board with built-in amplifier circuit?
  - Contacted Hamamatsu about this: they only offer modules with SiPM attached



Array board design by P. Bloser (U. of New Hampshire)


#### Readout of SiPM + PMTs in liquid nitrogen

- May reduce noise enough to see SER from SiPM
- Considering feeding LED flashes into setup via fiber-optic cable

# BONUS FEATURES...

### Calculating gain

Converting SER to gain:



W. Foreman


9

### The noise problem


LEMO cable that brings -1.5V square pulse to LED flasher passes close to PMT and SiPM bases

- Induces ringing in output
- Makes it harder to use "narrow flash" mode of LED w/inductor to do SER calibrations – though we get good results using "afterglow" method (w/ no inductor) as seen in previous slide (see "Light readout update (July 18 2013)" on DocDB for more info)






### The noise problem




Able to reduce ringing slightly using some readily available cable shielding\*

(Will try using steel mesh shielding and/or modifying box for different cable entry)



\* Shielding may or may not be aluminum foil


### SiPM circuit



# Suggested amp: Hamamatsu C4890 PIN photodiode amplifier



### Testing the SiPM



Even with lots of light, where 2-3 phe peaks should show up, we see nothing

From tests being done at PAB, we learned:

 expected SiPM signal on order ~0.1 mV, so it will be impossible to see it above the noise